

Invacare® Tracer® EX2 Wheelchair



FEATURES

- Durable, low-maintenance, triple chrome-plated carbon steel frame is long-lasting
- Dual axle positions allow variation of seat-to-floor heights
- Offset permanent arms available in 16", 18" and 20" widths
- 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Upholstery is dual-embossed
- Meets RESNA standards for manual wheelchairs in product performance

SPECIFICATIONS

Model Numbers	TREX20, TREX26, TREX28
HCPCS	K0001/K0002
Chair Dimensions	
Seat width	16", 18", 20"
Seat depth	16"
Back height	16"
Seat-to-floor height	17.5", 19.5"
Overall Dimensions	
Overall width	23.25", 24.25", 25.25", 26.25", 27.25", 28.25
Folded width	13"
Overall depth w/o riggings	31.88"
Overall height	34"
Weight w/o riggings	36 lb.
Weight capacity	250 lb.
Limited warranty (mo.)	60 months
Casters and Wheels	
Caster height	8"
Wheel height	24"
Shipping	
Box weight	43 lb.
Box dimensions (LxWxH)	32" x 13" x 38"

Note: All specifications are approximate

Invacare® Tracer® EX2 wheelchair sets the benchmark in manual wheelchairs. A true dual axle position allows for repositioning the 24" rear wheels and 8" casters for adult and hemi seat-to-floor heights. Offset permanent arms offer full seat width up to 20" with a permanent arm, hemi-spaced footrest and full range of accessories. Incorporating the design and technology of the Invacare 9000 series, the Tracer EX2 wheelchair has interchangeable components with the 9000 series chairs.



Ask For It!