

## Elevate Upright Walker

HCPCS:



The Drive Elevate Upright Walker brings users increased balance and stability while promoting better posture and providing a more natural line of sight than traditional mobility aids. Height-adjustable upholstered forearm pads keep arms supported and level for comfortable use, while large, ten-inch wheels make it perfect for indoor and outdoor use.

- Elevate Upright Walker is the stand-up mobility aid that's easy to take on the go.
- Padded armrests and upright handle adjust to 44"
- Durable and comfortable padded seat
- Simple one-handed, side-to-side folding
- Includes removable zipper storage bag with strap
- Dual front fork for optimal stability and
- Loop lock brakes for added safety and control
- Adjustable length dual handbrake system
- Large 10 multi-terrain front wheels for
- Height adjustable backrest
- \*Sit to stand assist handles currently on blue model only\*

PRODUCTS	Description
UPWKRBL10	UPRIGHT WALKER, ROLLATOR, 1EA, BLUE
UPWKRGY10	UPRIGHT WALKER, ROLLATOR, 1EA, GREY
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	23 lb
Material	Aluminum
Seat Dimensions	18" (W) x 10" (D)
Handle Height	39" - 44"
Weight Capacity	300 lb
Seat To Floor Height	21"
Casters	Front: 10   Rear: 8
CA Prop 65	no
Handle Extension Range	4" - 6"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]