

Viper Jr. Wheelchair

HCPCS: E1236,E1236



- Easier to propel and transport
- Frame style eliminates seat rail guides and allows for custom backs and accessories
- Back height adjustable from 18"20"
- Carbon steel black powder-coated frame double coated for an attractive, chip-proof, maintainable finish
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Composite Mag-style wheels with aluminum hand rims are lightweight and maintenance free
- Polyurethane rear tires and casters offer superior performance, smoother ride and are lightweight
- 8" Front casters have three height adjustments
- Nylon, padded armrests provide added user comfort
- Flip-back arms allow for easy transfers
- Composite steel lined front armrest receiver is angled so the arm will lock into place every time
- The front mounted release lever can be flipped up or down to unlock flip-up arm, depending on the patient's needs
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Dual axle provides easy adjustment of seat height
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite foot plates with heel loops are durable and lightweight
- Comes standard with flip-up anti-tippers
- Comes with push-to-lock wheel locks
- Available in 12" and 14" seat widths

PRODUCTS

Description

L412DDA-ELR	Viper,12"FlpBk.DskArm,ELR,1/cs
L412DDA-SF	Viper,12"FlpBk.DskArm,SF,1/cs
L414DDA-ELR	Viper,14"FlpBk.DeskArm,ELR,1/cs
L414DDA-SF	Viper,14"FlpBk.DskArm,SF,1/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Product Weight	30 lb 31 lb
Product Dimensions	36" (L) w/ riggings x 20" 22.5" (W) x 32" (H)
Seat Dimensions	12" 14" (W) x 12" 14" (D)
Weight Capacity	250 lbs
Seat To Floor Height	15.5" - 17.5"
Carton Shipping Weight	39 lbs or less
Folded Dimensions	36" (L) w/ riggings x 12" (W) x 32" (H)
Back Height	16" - 18"

SPECIFICATIONS

CA Prop 65

Description

[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.]

ACCESSORIES

Description

STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS2018N	Seat RailExt Kit,20x18,Blk.Nylon,1/ea, EA
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS807	Anti-Tips,Adj.w/Whls,Blk,1pr/bx,10/cv, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA