



171



199



Deluxe 3 Wheel Steel Rollator, 7.5" Casters

Loop Locks

Item #	Color	UOM
171	Chrome	1/cs

- Comes standard with basket, tray and pouch. (Figure A)
- Basket conveniently folds for easy storage. (Figure B)
- Ships in ultra-compact carton. Comes with frame folded and compact. Pull trombone-style frame to lock in place and secure with lock knob.
- Chrome-plated steel frame.
- Lightweight, solid 7.5" casters for indoor or outdoor use.
- Handles easily adjust in height with self-threading knob. (Figure C)
- Special loop lock made of internal aluminum casting operates easily and ensures safety. (Figure D)
- Easy, one-hand folding.

specifications

Handle Height	31" - 38"
Unit Dimensions	25" (W) x 24" (L)
Carton Shipping Weight	19 lbs
Weight Limit	300 lbs.
Warranty	Limited Lifetime

accessories

815B	Replacement Basket	1/ea
------	--------------------	------

Winnie Lite Supreme/Go Lite 3 Wheel Aluminum Rollator

7.5" Casters, Loop Locks

Item #	Color	UOM
199	Tan Plaid	1/cs

- Special loop lock made of internal aluminum casting operates easily and ensures safety. (Figure A)
- 7.5" casters with soft-grip tires are ideal for indoor and outdoor use.
- Comes standard with carry pouch.
- Handles easily adjust in height with self-threading knob. (Figure B)
- Easy, one-hand folding.

specifications

Handle Height	32" - 38"
Unit Dimensions	26" (W) x 22.75" (L)
Carton Shipping Weight	15 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime