Protekt Protopia™
"Ultra" Low Long Term Care Bed

THE **MOST RELIABLE AND COST-EFFECTIVE** BED IN ITS CLASS!

ULTRA LOW AND
HIGH HEIGHT RANGES





Extensive Travel Range 7"-31"
Expandable Length 80"-84"-88"

Lowest position of only 7"!
Helps mitigate injury from accidental roll outs

Bed Rail Options:


Quarter Rail


Pivot Assist Rail


QBar Assist Rail


Swing Rail

Color Options:


Mahogany


Amber Cherry


Oak

Features and Benefits:

- Heavy duty welded steel construction
- 8 function hand pendant
- Auto contour function
- Pair of quarter rails or assist bars and head & foot boards are included
- Integrated tool free length extension
- Bumper bar protects walls from damage
- Corner mattress retainers
- Quick release head and foot boards
- (8) 3" casters: 4 locking and 4 swivel
- **5 YEAR WARRANTY ON ELECTRICAL COMPONENTS**
- 15 Years on frame (1 Year on all other parts)

Model	Height		Length Expansion	Weight Capacity
	Min	Max		
Protopia™ #PMLTC36	7"	31"	80"-84"-88"	500 LBS.

Proactive medical products
270 Washington Street, Mount Vernon, NY 10553
Tel: 855-BE-PROACTIVE (855-237-7622) * **Fax:** 845-205-6005
www.proactivemedical.com * sales@proactivemedical.com

Protekt Protopia™
"Ultra" Low Long Term Care Bed

ITEM #	H/F Board Color	Rail Type
PMLTC36	No H/F Boards	No Rail
PMLTC36-MGY	Mahogany	No Rail
PMLTC36QR-MGY	Mahogany	Quarter Side Rail
PMLTC36QB-MGY	Mahogany	Qbar Assist Bar
PMLTC36PR-MGY	Mahogany	Pivot Assist Bar
PMLTC36SR-MGY	Mahogany	SoftTouch Swing Rail
PMLTC36-CHRY	Amber Cherry	No Rail
PMLTC36QR-CHRY	Amber Cherry	Quarter Side Rail
PMLTC36QB-CHRY	Amber Cherry	Qbar Assist Bar
PMLTC36PR-CHRY	Amber Cherry	Pivot Assist Bar
PMLTC36SR-CHRY	Amber Cherry	SoftTouch Swing Rail
PMLTC36-OAK	Oak	No Rail
PMLTC36QR-OAK	Oak	Quarter Side Rail
PMLTC36QB-OAK	Oak	Qbar Assist Bar
PMLTC36PR-OAK	Oak	Pivot Assist Bar
PMLTC36SR-OAK	Oak	SoftTouch Swing Rail

REPLACEMENT RAILS	
ITEM #	**Full Description**
PMLTCQRAIL	Quarter Side Rail f/ LTC Bed (pr)
PMLTCQBAR	Qbar Assist Bar f/ LTC Bed (pr)
PMLTCPRAIL	Pivot Assist Bar f/ LTC Bed (pr)
PMLTCSRAIL	SoftTouch Swing Rail f/ LTC Bed (pr)

REPLACEMENT HEAD & FOOT BOARDS	
ITEM #	**Full Description**
PMLTCHFB-MGY	Head & Foot Boards f/ Non-Expandable LTC Bed - Mahogany
PMLTCHFB-CHRY	Head & Foot Boards f/ Non-Expandable LTC Bed - Amber Cherry
PMLTCHFB-OAK	Head & Foot Boards f/ Non-Expandable LTC Bed - Light Oak

Proactive medical products

270 Washington Street, Mount Vernon, NY 10553
Tel: 855-BE-PROACTIVE (855-237-7622) * **Fax:** 845-205-6005
www.proactivemedical.com * sales@proactivemedical.com